
Earlier and Higher Fusion Rate

NOVOSIS

rhBMP-2 Recombinant Human
Bone Morphogenetic Protein 2

0.5mg | 1.0mg | 3.0mg

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
Ref.

1. MR Urist and BS Strates, Bone morphogenetic protein. *J Dent Res*, 1971. 50(6): p. 1392-406.
2. TK Sampath and AH Reddi, Dissociative extraction and reconstitution of extracellular matrix components involved in local bone differentiation. *Proc Natl Acad Sci U S A*, 1981. 78(12): p. 7599-603.
3. TK Sampath and N Muthukumar, et. al., Isolation of osteogenin, an extracellular matrix-associated, bone-inductive protein, by heparin affinity chromatography. *Proc Natl Acad Sci U S A*, 1987. 84(20): p. 7109-13.
4. JM Wozney and V Rosen, et. al., Novel regulators of bone formation: molecular clones and activities. *Science*, 1988. 242(4885): p. 1528-34.
5. KL Ong and ML Villarraga, et. al., Off-label use of bone morphogenetic proteins in the United States using administrative data. *Spine (Phila Pa 1976)*, 2010. 35(19): p. 1794-800.
6. WHO Drug Information, Vol. 28, No 1, 2014, List 71
7. NOVOSIS 허가증
8. M Wu and G Chen, et. al., TGF-beta and BMP signaling in osteoblast, skeletal development, and bone formation, homeostasis and disease. *Bone Res*, 2016. 4: p. 16009.

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

- * MTD : Maximum Tolerated Dosage
- * NOAEL : No-observed-adverse-effect level

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
Study type	Species/Strain	Route	Duration	Dose (mg/kg)	Relevant Findings
Single Dose Toxicity	Rat/Sprague-Dawley (male)	IV	single	0, 0.7, 1.8, 7	• MTD > 7 mg/kg
Repeated Dose toxicity	Rat/Sprague-Dawley (male & female)	IV	2 weeks	0, 0.05, 0.18, 0.5	• No observation of DWP431 related toxicity • NOAEL > 0.5 mg/kg

Ref.

1. JH Lee and SJ Jang et. al., Expression, Purification and Osteogenic Bioactivity of Recombinant Human BMP-2 Derived by Escherichia Coli. *Tissue Engineering and Regenerative Medicine*, 2011. 8(1): p. 8-15.
2. Data on file.
3. JH Lee and EN Lee, et. al., Acute Intravenous Injection Toxicity Study of Escherichia coli-Derived Recombinant Human Bone Morphogenetic Protein-2 in Rat. *Asian Spine J*, 2014. 8(2): p. 113-8.
4. JH Lee and EN Lee, et. al., The short-term effects of repetitive E. coli-derived rhBMP-2 administration through intravenous injection in rats. *Drug Chem Toxicol*, 2014. 37(1): p. 40-7.

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


Ref. 1. JH Lee and SJ Jang et. al., Expression, Purification and Osteogenic Bioactivity of Recombinant Human BMP-2 Derived by Escherichia Coli. *Tissue Engineering and Regenerative Medicine*, 2011. 8(1): p. 8-15.

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
*H&E staining : Haematoxylin and eosin staining

Bone Volume(mm³)

Ref. 1. Data on file.

NOVOSIS Carrier 탐색 1 Hydroxyapatite(HA)

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
*HAGE: Hydroxyapatite Granule Extralarge

Sample #	Macro porosity (%)	Interconnectivity (%)
1	71.6	99.9
2	70.2	99.5
3	65.9	99.6
Mean \pm SD	69.2 \pm 3.00	99.7 \pm 0.23

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

Ref. 1. JH Lee and SJ Jang et. al., Expression, Purification and Osteogenic Bioactivity of Recombinant Human BMP-2 Derived by Escherichia Coli. *Tissue Engineering and Regenerative Medicine*, 2011. 8(1): p. 8-15.
2. Data on file, BA06-DQ-REP-004, 합성재료이식용뼈 가공용 및 기공 연결성에 관한 보고서

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• Union rate(%)


Ref. 1. JH Lee, et al., Comparative study of fusion rate induced by different dosages of Escherichia coli-derived recombinant human bone morphogenetic protein-2 using hydroxyapatite carrier. *The Spine Journal*, 2012 Mar; 12(3):239-248.

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

• Union rate within cages(%)

• Overall union rate(%)


Ref. 1. CJ Hwang and JH Lee, et. al., Evaluation of the efficacy of Escherichia coli-derived recombinant human bone morphogenetic protein-2 in a mini-pig spinal anterior interbody fusion model. *Bone Joint J*, 2013. 95-B(2); p. 217-23.

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



Mineralized collagen fiber at bone tissue

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Ref.

1. WJ King and PH Krebsbach, Growth factor delivery: how surface interactions modulate release in vitro and in vivo. *Adv Drug Deliv Rev*, 2012. 64(12): p. 1239-56.
2. W Friess and H Uludag, et. al., Characterization of absorbable collagen sponges as rhBMP-2 carriers. *Int J Pharm*, 1999. 187(1): p. 91-9.
3. S Schaefer and R Detsch, et. al., How Degradation of Calcium Phosphate Bone Substitute Materials is influenced by Phase Composition and Porosity. *Advanced Engineering Materials*, 2011. 13(4): p. 342-350.
4. J Louis-Ugbo and HS Kim, et. al., Retention of 125I-labeled recombinant human bone morphogenetic protein-2 by biphasic calcium phosphate or a composite sponge in a rabbit posterolateral spine arthrodesis model. *J Orthop Res*, 2002. 20(5): p. 1050-9.
5. ML Bouxsein and TJ Turek, et. al., Recombinant human bone morphogenetic protein-2 accelerates healing in a rabbit ulnar osteotomy model. *J Bone Joint Surg Am*, 2001. 83-A(8): p. 1219-30.
6. W Friess and H Uludag, et. al., Characterization of absorbable collagen sponges as rhBMP-2 carriers. *Int J Pharm*, 1999. 187(1): p. 91-9.
7. PQ Ruhe and OC Boerman, et. al., In vivo release of rhBMP-2 loaded porous calcium phosphate cement pretreated with albumin. *J Mater Sci Mater Med*, 2006. 17(10): p. 919-27.
8. J Tazaki and M Murata, et. al., BMP-2 release and dose-response studies in hydroxyapatite and beta-tricalcium phosphate. *Biomed Mater Eng*, 2009. 19(2-3): p. 141-6.
9. EJ Carragee and EL Hurwitz, et. al., A critical review of recombinant human bone morphogenetic protein-2 trials in spinal surgery: emerging safety concerns and lessons learned. *Spine J*, 2011. 11(6): p. 471-91.
10. X Huang and CS Brazel, On the importance and mechanisms of burst release in matrix-controlled drug delivery systems. *J Control Release*, 2001. 73(2-3): p. 121-36.
11. X Dong and Q Wang, et. al., Understanding adsorption-desorption dynamics of BMP-2 on hydroxyapatite (001) surface. *Biophys J*, 2007. 93(3): p. 750-9.
12. J Patterson and R Siew, et. al., Hyaluronic acid hydrogels with controlled degradation properties for oriented bone regeneration. *Biomaterials*, 2010. 31(26): p. 6772-81.
13. T Boix and J Gomez-Morales, et. al., Adsorption of recombinant human bone morphogenetic protein rhBMP-2m onto hydroxyapatite. *J Inorg Biochem*, 2005. 99(5): p. 1043-50.
14. EJ Carragee and EL Hurwitz, et. al., A critical review of recombinant human bone morphogenetic protein-2 trials in spinal surgery: emerging safety concerns and lessons learned. *Spine J*, 2011. 11(6): p. 471-91.
15. JM Toth and SD Boden, et. al., Short-term osteoclastic activity induced by locally high concentrations of recombinant human bone morphogenetic protein-2 in a cancellous bone environment. *Spine (Phila Pa 1976)*, 2009. 34(6): p. 539-50.
16. JD Boerckel and YM Kolambkar, et. al., Effects of protein dose and delivery system on BMP-mediated bone regeneration. *Biomaterials*, 2011. 32(22): p. 5241-51.

The Clinical Trial for Evaluation of Efficacy and Safety of Novosis in Posterolateral Spinal Fusion¹

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