
## **Product Specifications**



				V I USII	er Length 187c
Туре		Catalogue No.	Primary Coil Diameter (inch)	Secondary Coil Diameter* (mm)	Primary Coil Length (cm)
		392-0304	0.010	3 -	4
	ExtraSoft (ES)	392-0306			6
		392-03H8		3.5	8
	Soft (S) MediumSoft (MS)	392-0408		4	8
		392-0410			10
		392-0510	0.012	5 -	10
		392-0515			15
		392-0620		6	20
		392-0730		7	30
Complex -		392-0840	0.014	8	40
		392-0940		9	40
		392-1040		10	40
		392-1250		12	50
	Standard	392-1450		14	50
	(STD)	392-1650		16	50
	SilkySoft (SS)	390-0101	- - - - - - - - - - - - - - - - - - -	1	1
		390-0102			2
		390-0103			3
		390-01H2		1.5	2
		390-01H3			3
$\sim$		390-0202		2	2
Complex		390-0203			3
		390-0204			4
		390-0206			6
		390-0208			8
		390-02H4		2.5	4
		390-02H6			6
		390-0304		3	4
		390-0306			6
		390-0308			8
Complex Infini	SilkySoft (SS)	393-020304	0.012	2-3	4
		393-020306			6
		393-020310			10
		393-030515		3-5 -	15
		393-030520			20

🗸 Pusher Length 187cm



 $\ast$  The secondary coil diameter means referenced diameter of a space to be embolized.

Туре	Catalogue No.	Description
EDG v4	347-404	ELECTRO DETACH GENERATOR v4

Do not re-use

Consult instructions for use **STERILEED** Sterilized using ethylene oxide

Wholesaler B. BRAUN KOREA

13F West Wing, Posco Center, 440 Teheran-ro, Gangnam-gu, 06194, Seoul, Korea Tel. +82.2.3459.7800 I Fax +82.2.3453.3060 www.bbraun.co.kr Manufacturer KANEKA CORPORATION 3-18, 2-Chome, Nakanoshima, Kita-ku, Osaka-city, OSAKA, 530-8288 JAPAN Tel. +81.6.6226.5256 https://www.kaneka-med.jp/english/



Aesculap Neurosurgery











## Framing that is stable & easily fits into aneurysms

• Designed to sit in open space

The distal coil shape spreads along the aneurysm wall & the proximal coil shape provides support from the inside.

• Coil designs optimized according to size

Combination of Primary Coil Diameter & length for easy delivery & stability.



### Complex SilkySoft



## Pursuit of softness, SilkySoft

### • 3D & Helical Hybrid

Straight, in-line structure coil design searches for space while entwined with other coils.

## • Softest coil in our product history

Soft coils & the Waving SR makes tight packing & minimal catheter kick-back possible.





# Complex Infini SilkySoft

• 12 Coil Type SilkySoft A soft coil with a 0.012" Primary Coil Diameter for volume maximization. • Filling coil that moves into space with ease of the space and pack them.



EDGv4 indicates by Sound and Light to signal the best Detach Point even when position of 2nd marker of Microcatheter is not detectable under X-Ray Fluoroscopy.



- Large loops that can fold tightly, configure according to the size
- Greater volume gain performance with a single coil
- Due to the Waving SR and softness of the SilkySoft, the coil can fold compactly enabling choice of longer coils possible.